

Empowering through Communication

Augmentative Communication Products Assistive Technology & Software

Behind the need to Communicate is the need to share.

Behind the need to share is the need to be Understood.

> Leo Rosten-American comic novelist (1908 - 1997)

TABLE OF CONTENTS

More Choices than ever in the Chat Family

Welcome to the 2012 Saltillo Products Catalog. I invite you to turn the page to see our NOVA chat devices, an exciting new line of devices based on the Android platform.

As the Chat hardware options have expanded, so have the vocabulary options. In fall 2011, we introduced ChatPower which is a word-based vocabulary that is used to create sentences 'on the fly.' ChatPower was developed by Nancy L. Inman, M.A.T., CCC-SLP specifically for the NOVA chat 7 & 10. ChatPower 24 & 42 were introduced for both NOVA chat 7 & 10 while ChatPower 80 is exclusively for NOVA chat 10. This past year we also introduced MultiChat Adult. There are multiple forms of communication a vailable in this program; sentences, phrases, individual words, recordings for story-telling and visual scenes.


In addition to product enhancements, much of our focus is still spent on providing the best service and support resources for our customers. In 2011 we launched a new website (Saltillo.com) which includes an easy to navigate and very comprehensive knowledge base. Our network of regional consultants continues to grow which allows us to provide personal training and support. We've also became very active in the social media arena with Saltillo Facebook, Twitter, and YouTube groups.

If you wish to discuss any of the products in the catalog or need guidance in making a selection, please call us at 1-800-382-8622 or email us at support@saltillo.com. We will gladly assist you with any questions that you have.

Sincerely,

Dave Hershberger Vice President









The NOVA chat Family Living Life to its Fullest.

NOVA chat 7 Chat Software with a 7" Display



NOVA chat 7



NOVA chat 10 Chat Software with a 10" Display



NOVA chat 10

NOVA chat 7

Living Life to its Fullest.

IVONA synthesized speech (including teen and child voice) or optional Acapela speech

SymbolStix or optional PCS symbols



Size: 7.48" x 4.74" x 1.62" Weight: 1.44 lbs.

NOVA chat 7

4 Interchangeable color shells

NOVA chat 7 is a communication device that offers the well known Chat software on an Android platform. Like all "Chats," NOVA chat 7 offers a durable yet sleek portable design, with features that include a 7" display, switch scanning, Social Chat, new IVONA speech synthesizer, and many other Chat features.

NOVA chat 7 includes thousands of SymbolStix™ symbols, PCS™ symbols are available as an option. ChatPower 24 & 42 as well as MulitChat 15 are a couple of the vocabulary programs included in NOVA chat 7.





Vocabulary files for both literate and emerging literate clients

Large 7.0 " WSVGA capacitive touch screen

Configurable orientation (portrait or landscape)

Built-in stand





Chat Software

Chat software is a feature-rich application that provides the power behind the Chat products. This software contains the functionality and flexibility to meet the unique communication needs of individuals with varying needs and preferences. The included vocabulary files are an excellent starting point for many children and adults who need a portable communication device.

Pages

Like many dynamic screen communication devices, Chat software allows pages to be linked together to form a communication environment. Pages may be organized by activities, as you will find in the myQuickChat vocabularies, or organized in a way that leads you seamlessly through the construction of a sentence, as in ChatPower.

You have the flexibility to modify each page to meet the visual, access and language needs of each individual. You can select the number of buttons per page, the appearance of each button and the action that occurs when a button is pressed.

Vocabulary Files

All Chat devices come with several vocabulary files. These page sets are often great launching pads for individuals just beginning to use their communication device. The various vocabulary files are designed to meet a broad variety of communication needs. See pages 8-10 for details.

Symbols and Images

Symbols are often used on buttons to represent words and sentences, particularly for individuals who have difficulty with reading. Chat software comes with pre-loaded symbols for personalizing pages. SymbolStix is the default symbol set that comes with the Chat software. PCS symbols are an optional alternative.

Sometimes it is better to use a photo on a button instead of a symbol. For example, you may want to use photos of your classmates or family members. NOVA chat's offer built in cameras while computer programming software allows you to use any photo or image that you have on your computer as a picture for a button. Chat software also gives you the opportunity to use a photograph as the background for a page. On these pages, you can create buttons or hotspots to allow a beginning communicator to choose vocabulary from the context of a visual scene.

Spelling and Word Prediction

All Chat products include a variety of keyboards for spelling. Typed words can be spoken by the internal voice synthesizer. All 'Plus' and 'Dedicated Plus' models also include word prediction and abbreviation expansion.

Speech

All Chat devices include both digitized and synthesized speech. Digitized speech is recorded by a sibling, friend, etc. Synthesized speech is computer generated 'text-to-speech.' For the specific speech options available with each device, see page 16.

Windows Editing Software

You can easily customize and add vocabulary directly on the Chat device and/or use the included Editor software, which allows editing and customizing on a Windows computer. The Editor also provides an easy means of importing photos and other graphics to be used as symbols on the Chat device.

Running other Applications

Since Chat products are built on computers, other software applications can run on these hardware platforms. Chat software allows applications to be launched from a communication page button. Note: This feature is prohibited in dedicated models by Federal and state funding regulations.

Scanning

In addition to touch screen access, NOVA chat has a full complement of visual and auditory scanning capabilities. The user interface allows easy adjustment of various timing aspects, scanning patterns and switch input options to best take advantage of each individual's unique capabilities.











Chat Vocabulary

Chat devices ship with a number of default vocabulary files, or page sets. The vocabulary and organization of each vocabulary file was designed to address specific communication needs. All vocabulary files are fully customizable and can be personalized for each individual. The descriptions on the next three pages highlight a few of the available vocabulary files. For specific vocabulary included with each device, see page 16.

ChatPower

ChatPower is a word-based vocabulary using 24, 42 or 80 buttons per page. ChatPower was developed by Nancy L. Inman, M.A.T., CCC-SLP specifically for the NOVA chat.

ChatPower uses high-frequency core words to facilitate quick and easy sentence generation. It also includes nouns and adjectives that are organized in logical categories. Grammatical word classes have been included. As a sentence is being built, you often find that the next word you want to say can be spoken with one or two button presses. A spelling/word prediction page is used to spell words that are not included in the vocabulary set.

Being a word-based vocabulary, ChatPower can be used by a wide age group. Nouns, adjectives and verbs can be added and customized to meet the particular needs of the individual. Literacy skills are necessary to use ChatPower, as picture support has only been provided for picture producing words.

ChatPower24

ChatPower24 is a word-based vocabulary that allows for easy intuitive communication. The program has been designed specifically for this portable and readily available platform. It is a generative language system that will be familiar to those using WordPower on other systems, and easy to learn for new users.



(D) what I am thin

ChatPower24 with Phrases

ChatPower24 with Phrases is also a word-based vocabulary program, but phrases and sentences starters have been incorporated in order to make language generation easier for those that may have difficulty with a word-based program. Core language and category-based words are included in this program.

ChatPower42

ChatPower42 is a word-based vocabulary that allows for easy and intuitive communication. ChatPower42 is a generative language system that will be familiar to those using WordPower on other systems, and easy to learn for new users.



ChatPower42 Basic

ChatPower42 Basic is a simplified version of ChatPower42. This word-based system remains rich in core vocabulary, and it may be appropriate for beginning users or those that need a more basic vocabulary design. It is designed to make communication fast and easy.

ChatPower80

ChatPower80 created specifically for the NOVA chat 10 combines spelling and word prediction with core vocabulary all at once.

AACelerate with myVoice

AACcelerate with myVoice is a unique suite of AAC vocabulary sets for the Chat family targeting the communication needs of a broad spectrum of individuals. This collection, which is intuitively woven into several different core vocabulary sets and phrase-based sets, incorporates the AAC vocabulary research of Dr. David Beukelman, Dr. Diane Bryen and Barbara Collier.

myCore

myCore is a combination of core and phrase based vocabulary geared toward individuals with literacy skills. myCore is only available for the ALT-Chat.

myQuickChat

myQuickChat is an introductory communication system geared toward child and adult AAC users with complex communication needs. It offers a high frequency, phrase-based communication system in an easy-to-use, yet very engaging format. The progressive system offers a variety of topics for everyday needs and conversation and provides immediate and more successful communication exchanges within a variety of settings. myQuickChat is available for both children and adults in 4, 8 and 12 locations, each with an identical 16 location Support Master Home Page.

VocabPC™

VocabPC was designed by Gail Van Tatenhove, PA, MS, CCC-SLP, for adults and adolescents with developmental disabilities. Vocabulary in VocabPC is arranged as carrier phrases, interactive sentences, activity vocabulary and naming words. VocabPC uses a 12 location page layout. A tour guide for VocabPC is included on the NOVA chat Editor CD. This tour guide describes the users the author was working with and provides rationales and teaching support for VocabPC.

Let's pla	y again.	
8-8	2-8	R.
8	222	8-8
-	27	- TO
業	1	



- P	w	ت	<u> </u>	<u> </u>	Y	-		0	P
. 1	•		đ						
	110		1	<.	¥.	ь		m	
1	-	10	n	T	-	11	TT.	12	
~	6	34-	100	4	-	+	ant		1
		44	12	127	stat "	Aur .	-	4	
-	-	And and Post	57	(F)	1816	1000	-	with .	
It's	VOLU	r turi	2			_			
	100	Curr	**	-	-	-	-		-
	120	100			ALC: N	22		- 1	2

			sv t	-		Canal -	ante B			40		
-	Ant	-	y99		-	-	-		nete:	-	141	1
1	-		hed.	-	pet.	make	248	-	adated.		Ante	
-	-		Sec.		-	ment		-		-	-	met.
-	-	ant.		call		phage.	Farm	and.	-	for.	-	Date
.be	bim.	enald	and the	-	MA	aut.	-		100	hum		then.
aba	feet.	-04	-	-	-		-		het.	-	1	
Day.	tien		-P		-		waters	Base	54-	See.	-	-
	-	have	mould	feet,	-	Jahn.	-	-	1	-	-	

Hey there.			
QR	-	92-1	
	-	?	
?	0.0		

Chat Vocabulary cont.

MultiChat 15

Two versions of this 15 location vocabulary were created, one for school-aged individuals with emerging language skills and one for adolescents or adults. There are multiple forms of communication available in this program; sentences, phrases, individual words, recordings for story-telling and visual scenes. Also included are interactive play and reading pages, social pages and a News-2-You starter page for the weekly newspaper's vocabulary. There are symbols on every button with the exclusion of several core words.

...

4-Basic

Just like it sounds, a basic vocabulary option with 4 buttons per page.

Essence

Essence is specifically designed for adults with acquired speech disorders. The home page contains a spelling keyboard in addition to quick links to pages of commonly used phrases for communication with those you regularly interact with in your family and community. Other pages contain commonly used vocabulary that may be needed for a visit to the doctor, using public transportation, going out to eat, etc. Each page is completely customizable to meet the specific vocabulary needs of each individual.

Spelling

Spelling is a keyboard page set with four word predication buttons and a few pre-stored phrases.

Editor Software

NOVA chat Editor is the Chat software that is able to run on a Windows desktop or laptop computer. Although it can be done on the device, the Editor provides family and/or professionals a way to customize vocabulary for the device user when it is convenient.

Note: Having the vocabulary on a second computer is always recommended as a backup.



0

67

Symbol Options for NOVA chat and ChatPC devices

Saltillo Corporation offers two options for symbols in Chat products.

SymbolStix

The default symbol set for NOVA chat, ChatPC-Silk, and ALT-Chat devices is SymbolStix.



SymbolStix is a dynamic symbol set created by N2Y for use in its Internet weekly newspaper, News-2-You, and its online curriculum, Unique Learning System. Generic activities and persons are depicted as lively, vibrant stick Figures - stick figures with "an attitude." (Actual persons incorporated in SymbolStix are drawn as easily recognizable characterizations.) The current suite includes over 4,000 individual drawings/icons.





Examples of SymbolStix symbols.

Note: SymbolStix is not available in ChatPC-4

For more information on SymbolStix, visit www.symbolstix.com.


PCS


You can choose PCS symbols in additioin to SymbolStix when you order a Chat device. PCS symbols are very popular in classrooms and clinics throughout the world. These symbols have been created by Mayer-Johnson and are used in Mayer-Johnson's BoardMaker and Speaking Dynamically applications. Several thousand symbols represent vocabulary for a wide variety of school, home and work activities. PCS symbols can be selected for Chat devices as a \$375 option.





Examples of PCS symbols.







The ChatPC Family

ChatPC-Silk Top performance in a sleek package

ChatPC-4 Rugged and Versatile



0 like it better than the first 5ame 123

ChatPC-Silk

ALT-Chat ChatPC functionality on a 7" display





ChatPC features.

ALT-Chat includes thousands of SymbolStix symbols, PCS symbols are available as an option. Installed vocabulary programs include PalmChat, myVoice, VocabPC, 16 Primary and 16 Basic.





ChatPC-Silk

The sleekest and most powerful ChatPC ever!

0

an tell you what I am thinking

NeoSpeech, Loquendo or Acapela speech synthesis

Accessed by finger touch

Holds several hours of recorded speech



Can hold numerous vocabulary files, each containing hundreds of pages

> Large 4.0" (diagonal) screen

Removable protective rubber housing

> Software and hardware volume

> > Long-life rechargeable Li-lon battery

Size: with protective housing: 1.6" x 3.25" x 5.6" without protective housing: 1.4" x 2.95" x 5.25"

> Weight: with protective housing: 12.5 oz. without protective housing: 9 oz.

> > 12 to 36 month warranty

ChatPC-Silk

The ChatPC-Silk represents the ultimate in portability without sacrificing capability. Although it boasts the largest screen for any PDA-based communication device, it only weighs 12.5 oz with the protective rubber housing and 9 oz without. The ChatPC-Silk can quickly be inserted into the rubber housing to provide extra protection, or removed to reduce size and weight.

> ChatPC-Silk includes thousands of SymbolStix symbols, PCS symbols are available as an option. Installed vocabulary programs include PalmChat, VocabPC, 16 Primary and 16 Basic.

ChatPC-4[™]

The most rugged member of the ChatPC family!





A Comparison of the Hardware Platforms

All Chat products share many features and capabilities but there are also aspects that differentiate them. The chart below will help you decide which device best meets your needs.

	NOVA chat 7	NOVA chat 10	ALT-Chat	ChatPC-Silk
Strengths	Portability and power in a sleek durable package.	The newest NOVA chat member combines the powerful Chat Software on a 10" device.	Its large screen and scanning features make ALT-Chat a good option for those who cannot easily access the smaller Chat device.	Sleek, superior sound quality and all round, the most portable member of the Chat products.
Size	7.5" x 4.7" x 1.6"	10.5" x 7.5" x 1.5"	7.3" x 7.3" x 1.8"	3.3" x 5.6" x 1.6"
Weight	1.5 lbs	2.2 lbs	2.5 lbs	12.5 oz
Screen Size and Type	7" WSVGA capacitive touch	10.1" WSVGA capacitive touch	7" WVGA	4" VGA
Available Speech Synthesizers	IVONA Acapela (optional)	IVONA Acapela (optional)	NeoSpeech Acapela	NeoSpeech Loquendo Acapela
Symbol Sets	SymbolStix, PCS (optional)	SymbolStix, PCS (optional)	SymbolStix, PCS (optional)	SymbolStix, PCS (optional)
Included Vocabulary Sets	4 Basic ChatPower 24 ChatPower 24 with Phrases ChatPower 42 ChatPower 42 Basic Essence MultiChat 15 Adult MultiChat 15 Student myQuickChat Spelling VocabPC	4 Basic ChatPower 24 ChatPower 24 with Phrases ChatPower 42 ChatPower 42 Basic ChatPower 80 Essence MultiChat 15 Adult MultiChat 15 Student myQuickChat Spelling VocabPC	PalmChat VocabPC Primary 3x6 Saltillo Sampler Spelling Essence MultiChat 15 Student myQuickChat myCore	PalmChat VocabPC 16 Basic 16 Primary Saltillo Sampler Spelling MultiChat 15 Student myQuickChat
Available Languages	English	English	English-Spanish English-German	English-Spanish English-German
Models and Prices	Standard \$3495 Plus \$3995 Dedicated+ \$4995	Standard \$3895 Plus \$4395 Dedicated+ \$5395	Standard \$3495 Plus \$3995 Dedicated+ \$4995	Standard \$2395 Plus \$2995 Dedicated+ \$3995

Ordering and Price Information

When ordering a Chat device there are a few items you will need to consider so that your client or family member receive the most appropriate functionality to meet their needs. Here are the steps to make that happen.

STEP 1:

Choose a model Chat devices are available in three

models. Models differ from each other in features and capabilities as well as in the warranty period. Below is a description of each model.

Standard Model

The standard models are designed for individuals who will primarily be using symbols for communication. Although these models have keyboard pages for spelling, they do not include word prediction or abbreviation-expansion. Standard models come with a one year warranty. Standard models are "open" systems, meaning you can also run other platform specific applications on these devices.

Plus Model

The Plus models are designed for individuals with existing or emerging literacy skills. In addition to all of the features of the Standard models, Plus models come with 10,000, 30,000 and 60,000 word dictionaries. Numerous keyboard pages with and without spelling prediction are included in the pre-stored vocabulary files. Plus models also include abbreviationexpansion, a feature that allows you to encode a word, phrase or sentence under a short sequence of letters. Plus models also provide the capability to load and save text files. Plus models are "open", meaning you can also run other platform specific applications on these devices. Plus models come with a one year warranty.

Dedicated Plus Model

The guidelines for Medicare and for Medicaid in many states require that devices purchased with funds from these organizations are "dedicated" devices. Dedicated devices are devices that can only be used as communication devices, i.e. no other applications can be run on the devices. The Dedicated Plus Chat Models meet the funding requirements of these organizations. Dedicated Plus models include all of the features of the Plus models and come with a three year warranty. Dedicated Plus models also include a carrying case.

STEP 2: Choose a speech synthesizer

When ordering a NOVA chat 7 or NOVA chat 10 the IVONA speech synthesizer is included. As an option you can add the Acapela speech synthesizer, the voices will be added to existing IVONA voice options. The cost to include Acapala voices to NOVA chat is \$95. Below is a description of the NOVA chat synthesizers.

- > IVONA (IVONA) is a leader in the field with natural voice quality and inflection. IVONA voices include a child's voice, teenage voice, 2 adult female US English voices, one adult male US English voice, two adult female British voices and an adult male British English voice option.
- > Acapela (ACAP) voices include Graham: UK Male, Peter: UK Male, Ryan: US Male, Nizareng UK Male, Lucy: UK Female, Heather: US Female, Laura: US Female, Rachel: UK Female, Tracy: US Female, Kenny: US Boy (child), Nelly: US Girl (child).

When ordering a ChatPC or ALT-Chat, you need to specify a speech synthesizer. Below is a description of each synthesizer and its device availability.

- > Loquendo (LOQ) has 2 adult voices in US English, UK English, French, German, Italian and Spanish. Loquendo is an available option for ChatPC-Silk.
- > NeoSpeech (NEO) is the most natural sounding voice. NeoSpeech has 2 adult voices (male and female) in US English. NeoSpeech is an available option for ChatPC-Silk and ALT-Chat.
- > Acapela (ACAP) is available in English-Spanish or English-German variations. English voices include both adult and child voices. Acapela is an available option for ChatPC-Silk and ALT-Chat.

STEP 3: Specify Symbol Sets

By default, all Chat devices include SymbolStix symbols. PCS symbols can be specified by adding PCS to the end of the device model number and adding \$375 to the device price. For example ALT-NEO-BLU-PCS specifies PCS symbols for an ALT-Chat. See page 11 for a description of both symbol sets.

STEP 4: **Specify Bilingual devices**

NOVA chat 7 & 10 are only available in English, By default ChatPC-Silk and ALT-Chat have bilingual English-Spanish capabilities. Please note: Acapela or Loquendo speech synthesizers should be specified for devices that will be used as bilingual devices.

STEP 5: Choose color

NOVA chat devices come with 4 interchangeable color shells, the unit will ship with Gray unless another color is specified on the order. ChatPC-Silk has no color option. ALT-Chat is the only Chat device you need to specify a color, ALT-Chat is available in black, blue, and red.









Ordering and Price Information cont.

Below is an additional breakdown of ordering and pricing information for Chat devices. Refer to pages 16 and 17 for more detail on specific features. For clarity, we have color coded the part numbers in the description below. The hardware version is specified by the black font. The model is specified in by the blue font. Voice synthesizer is specified in green and color is specified in purple. Optional symbol and language indicators are show in red.

A 1998

11 12

STEP 2: Choose a speech synthesizer

ALT-NEO-BLA.....\$3495

- Use ACAP to specify alternate

- Use PCS to specify alternate

STEP 4: Specify Bilingual devices



NOVA chat 7 STEP 1: Choose a model STEP 2: Choose a speech synthesizer STEP 3: Specify Symbol Set NOVA chat 7, Standard model

- NOVA7-IVONA \$3495 - Use ACAP to specify optional synthesizer \$95 - Use PCS to specify optional
- symbol set.....\$375 NOVA chat 7, Plus model
- NOVA7-Plus-IVONA......\$3995 - Use ACAP to specify optional
- symbol set.....\$375
- NOVA chat 7, Dedicated Plus model NOVA7-DPlus-IVONA.....\$4995 - Use ACAP to specify optional synthesizer\$95
- Use PCS to specify optional symbol set.....\$375





STEP 2: Choose a speech synthesizer STEP 3: Specify Symbol Set NOVA chat 10, Standard model NOVA10-IVONA \$3895 - Use ACAP to specify optional synthesizer \$95 - Use PCS to specify optional symbol set.....\$375 NOVA chat 10. Plus model NOVA10-Plus-IVONA......\$4395 - Use ACAP to specify optional synthesizer \$95 - Use PCS to specify optional

STEP 1: Choose a model

symbol set.....\$375 NOVA chat 10. Dedicated Plus model NOVA10-DPlus-IVONA \$5395 - Use ACAP to specify optional synthesizer \$95 - Use PCS to specify optional

symbol set.....\$375

- Use RED to specify Red ALT-Chat, Dedicated Plus model ALT-DPlus-NEO-BLA \$4995 - Use ACAP to specify alternate
 - svnthesizer

ALT-Chat

STEP 1: Choose a model

STEP 5: Choose color

synthesizer

ALT-Chat. Plus model

synthesizer

STEP 3: Specify Symbol Set

ALT-Chat, Standard model

- Use PCS to specify alternate symbol set.....\$375

symbol set.....\$375

- Use BLU to specify Light Blue

- Use BLU to specify Light Blue

See page 33 for carrying cases,

accessories.

straps, mounting plates and other

- Use RED to specify Red

ChatPC-4 STEP 1: Choose a model STEP 2: Choose a speech synthesizer STEP 3: Choose color ChatPC-4, Standard model - Use LOO or DEC to specify alternate synthesizers - Use BLU to specify Blue case ChatPC-4, Plus model - Use LOQ or DEC to specify alternate synthesizers. - Use BLU to specify Blue case ChatPC-4, Dedicated Plus model CPC4-DPlus-NEO-BLA \$4195 - Use LOQ or DEC to specify


alternate synthesizers



- ChatPC-Silk STEP 1: Choose a model STEP 2: Choose a speech synthesizer STEP 3: Specify Symbol Set STEP 4: Specify Bilingual devices ChatPC-Silk, Standard model
- CPCSilk-NEO\$2395 - Use LOQ or ACAP to specify alternate synthesizers
- Use PCS to specify alternate symbol set.....\$375
- ChatPC-Silk, Plus model
- CPCSilk-Plus-NEO..... \$2995 - Use LOQ or ACAP to specify alternate synthesizers
- Use PCS to specify alternate
- symbol set.....\$375 ChatPC-Silk, Dedicated Plus model
- CPCSilk-DPlus-NEO \$3995
- Use LOQ or ACAP to specify
- alternate synthesizers - Use PCS to specify alternate
 - symbol set.....\$375



CPC4-NEO-BLA.....\$2495 CPC4-Plus-NEO-BLA..... \$3195




WordPower[™] Software by Inman Innovations, Inc. Simple choices. Fast responses.

WordPower software is a unique, word-based vocabulary design for an AAC device. WordPower combines the features of core vocabulary, spelling, and word prediction. The system takes advantage of the fact that a core vocabulary of just 100 words accounts for approximately 50 percent of words spoken. The core words of WordPower are categorized, colorcoded and alphabetized for easy access. WordPower uses these 100+ words in conjunction with spelling and word prediction. Picture WordPower is also a word-based, generative vocabulary for individuals needing picture support and/or larger cells for access. Picture WordPower consists of single-hit words, category-based words, and spelling with word prediction. Picture WordPower also now includes a Phrases and Categories version - a generative vocabulary option for individuals not able to use a strictly word-based system.

WordPower is available as software that runs on a variety of communication devices from Prentke Romich Company. Please note that the communication devices listed below must be purchased separately from Prentke Romich Company or one of their distributors.

WordPower[™] for ECO2

WordPower and Picture WordPower for the ECO are upgraded vocabularies that have been programmed to take advantage of the newest hardware and software innovations from Prentke Romich. The ECO2 is an integrated AAC device. WordPower and Picture WordPower are sophisticated vocabularies that have been designed to make speech generation easy and natural. In addition to these logical and intuitive vocabularies, WordPower and Picture WordPower now serve as the interface for the integrated ECO computer. The system operator can move smoothly between communicating with the device, navigating the internet, playing games on the computer, sending e-mail, and generating documents using WordPower vocabularies.



WordPower for ECO2

#WP-ECO\$350

WordPower[™] for Vantage Lite

WordPower and Picture WordPower vocabularies have been updated and revised for Prentke Romich's Vantage Lite. The vocabularies are now available in 45, 60 and 84-location configurations. There is a new 84-location version of Scanning WordPower.

WordPower for Vantage Lite

#3206\$350



WordPower[™] for Vantage Lite

WordPower[™] for ECO2



symbol set.....\$375 - Use BLU to specify Light Blue - Use RED to specify Red ALT-Plus-NEO-BLA \$3995 - Use ACAP to specify alternate Use PCS to specify alternate

WordPower[™] OnBoard

For those times when you need a manual communication board, WordPower OnBoard is the perfect choice. Based on the original WordPower vocabulary design by Nancy Inman, MAT, CCC-SLP, this sturdy, two-sided 9.5" x 11.5" board allows literate individuals to point to words with their fingers or a light pointer. The design allows it to be customized using blank cells to add words, phrases or sentences.

Featuring 100 of the most commonly used words, WordPower Onboard can also be used as back-up for those using another communication device. > Size: 11.5" x 9.5"

WordPower OnBoard #3202\$25



WordPower[™] OnBoard



SpringBoard Life[™] by Prentke Romich Company

Springboard Lite is a dynamic-screen digitized-speech communication device from Prentke Romich Company (PRC) that uses the same hardware as ALT-Chat. It has been updated with a widescreen format 7" touchscreen. Springboard Lite is a digitized device; the included vocabulary files have pre-recorded male and female voices. New recordings can easily be made when creating additional messages.

Vocabulary-Building Features

- > Five integrated versions of pre-stored vocabulary using 4-, 8-, 15-, and 32-locations with core vocabulary.
- > Pre-stored Activity Row categories
- > Six User Areas.
- > Visual Scenes for concept and language learning.

Easy Access Methods

- > 4-, 8-, 15-, and 32-location displays
- > Direct selection
- > Single- and dual-switch scanning
- > Joystick

Advanced Features

- > New "widescreen TFT display" offers improved picture quality.
- > Infrared environmental control functions to control common IR devices, such as lights and TV/ DVD players.
- > USB flash memory port and Secure Digital card slot enable faster, easier memory backup and transfer options.
- > Pre-loaded male and female recorded voices.

Customization Aids

- > Access to over 3,000 Minspeak® icons in the device for customizing vocabulary.
- > Over 8,500 optional PCS symbols for core vocabulary, Activity Rows, and Pages.
- > Ability to add your own graphics for symbols or visual displays

Device Specifications

- > Weight: 2 lbs. 8 oz.
- > Dimensions: 7.3"W x 7.3"H x 1.8"D
- > Display size: 7" WVGA (diagonal measurement)
- > Speech memory: Over 60 minutes of digitized speech capacity
- > Battery life: 6-8 hours with normal operation

SpringBoard Lite[™] and Bilingual SpringBoard Lite™

SpringBoard Lite English

- Includes:
- > Battery charger > Complete operator manual on
- CD-ROM
- > "Quick-Start" operator's guide
- > Unity[®] vocabulary guide
- > Starter curriculum
- > Two-year PRC warranty

Specify Color

(Black, Seafoam, Orange, Grape) #SBL-AEN\$2595

SpringBoard Lite Bilingual

Spanish-speaking users can communicate in either language at home, work, and school, with the device fully supported in Spanish the instruction manual is provided in both languages, and bilingual customer service support is available. Includes:

- > Battery charger
- > Complete English and Spanish
- operator manual on CD-ROM
- > "Quick-Start" operator's guide
- > Unity[®] vocabulary guide
- > Starter curriculum
 - > Two-year PRC warranty

Specify Color

(Black, Seafoam, Orange, Grape) #SBL-ASP\$2845





OPTIONS:



(Must be specified at time of purchase)



Note: Although Springboard uses the same hardware platform as the ALT-Chat, the software is very different. ALT-Chat has both digitized and synthesized speech. Springboard has only digitized (recorded speech). The included vocabulary is also different. Springboard vocabulary is based on Prenke Romich Company's Unity program. ALT-Chat is described on page 13.

I want to play





SpeakOut by permobil

SpeakOut is a portable communication device with a user-friendly QWERTY keyboard. The primary LCD display is used to compose messages, show word prediction options and navigate menus. The composed messages also appear on the rear screen, ideal for face-to-face communication.

SpeakOut uses Acapela synthesized speech. It includes both an adult male and adult female voice. A 10,000 word 'learning' word predictor can be used to assist with typing speed and spelling accuracy. Frequently used words or phrases can be stored under a key combination. You can also create documents that can be transferred to a computer or printed directly on an HP printer.

SpeakOut has a number of built-in tools that can be helpful throughout the day. Reminders can be set for important events or appointments, phone numbers can be stored in the included phone book and SpeakOut also includes a built-in calculator.



SpeakOut



Features

> Large, easy-to-read display

- > Acapela voices
- > Simple menus
- > QWERTY keyboard
- > 2 font sizes
- > Reminders
- > Phrase system
- > Word prediction with 10,000 words
- > Phone book
- > Large storage capacity
- > Text files can be transferred to a computer
- > Calculator
- > Time and date > Direct print to many Hewlett Packard
- printers
- > Amplified speaker
- > Includes canvas carrying case and shoulder strap

Technical Details

- > Rear LCD-Display, 40 x 2 characters
- > Front LCD-Display, 240 x 64 pixels
- > 1 USB for transferring information
- to PC > 1 USB for external devices, e.g.
- printer
- > Serial plug for IR-transmitter
- > Size: 9.7" x 7.7" x 2.0"
- > Operating time: 8 hours
- > Weight: 2.2 lbs.
- > Battery: Li-ion
- SpeakOut #SPKO-01\$3995



Scanning SpeakOut

are included.

Scanning SpeakOut



SpeakOut


The Scanning SpeakOut is a good option for those who cannot access a keyboard, but are able to activate a switch for scanning. The Scanning SpeakOut's keys cannot be pressed manually; it can only be operated via scanning or with the included USB keyboard. The USB keyboard is also required for certain setup features. A USB keyboard and USB hub

#SPKO-01SC\$4495

Accessories Keyguard - 3mm:

Isolates the keys, helping to stabilize
and position fingers or pointers. The
keyguard is made of high-quality
acrylic. 3mm (0.12") thick.
#SPKO-KG3\$235
Keyguard - 5mm:
Isolates the keys, helping to stabilize
and position fingers or pointers. The
keyguard is made of high-quality
acrylic. 5mm (0.20") thick
#SPKO-KG5\$235
Water Resistant Cover:
A clear, plastic rain cover that covers
the entire device.
#SPKO-RC \$115
Keyboard Protector:
A clear, plastic protector that covers
the keyboard, but not the display.
#SPKO-KP\$115
Hard Carrying Case:
Aluminum suitcase with foam
customized for the SpeakOut.
Protects the unit during transport.
Includes shoulder strap and keys for
locking the suitcase.
#SPKO-HCC\$165

Hard Carrying Case



ChatBox40[®] and ChatBox40–XT[®] Grow and learn with this versatile friend

· O M O

Is there a child or adult in your life or under your care that longs to communicate but cannot due to an inability to speak?

With the ChatBox40 and ChatBox40-XT, Saltillo offers two innovative products that open a world of communication. Using these durable and attractive products, individuals who feel the frustration of not being able to express themselves can find greater fulfillment in their lives and in their interpersonal relationships.

Using the ChatBox40, the press of a button (or one of various switch scanning options) results in communication! As an aid for literacy development, the spoken message can also appear on the text display. A single word, a complete sentence, or a recorded song can be readily expressed with one button activation.

Conversations can be held using either the forty key standard keyboard or the alternate keyboard which features ten large keys. Use the up and down arrows to access different levels where additional messages are stored.

Cat

For individuals with a visual impairment, the ChatBox40-XT offers auditory prompts.

Another great plus of the ChatBox 40-XT, which adds to its growth potential and flexibility, is a basic voice synthesizer. Unique messages can be created using either the pre-stored messages or by typing on the alphabet keyboard. The ChatBox40-XT may also be connected to a computer to imitate a regular keyboard.

With models available in eyecatching colors, a young child can grow with the ChatBox40 to reach a level of basic literacy, while older individuals can expand their abilities. Both devices are highly adaptable, and their longevity and durability make them very cost-effective choices.

ChatBox40 (Gray) #2600GRAY\$895

Key ChatBox40 Features

Digitized Speech – Thirty minutes of high-quality, recorded speech can be used as needed across the buttons in the 10 overlays. A single key can be programmed with almost two minutes of recorded speech. Lights (LEDs) on each key provide the user with a visual cue of stored messages.

Text Generation – Digitized messages can be programmed to include text to appear in the 2-line LCD Display. Using this feature, an individual can assemble a message using multiple keys (e.g., "I want ... to eat"). The message can then be spoken as a single message by selecting the Speak Message button. On the ChatBox40-XT, the text message can also be spoken with the built-in speech synthesizer.

Scanning Options - The ChatBox40 supports single switch, dual switch, and joystick scanning. Scanning options include row/column, linear, directional, latched, direct, and step scanning. All scanning modes are visually supported with lights on the keys. The scanning time is variable from .1 - 10 seconds. Auditory prompting (row and key) is available on the ChatBox40-XT.

Various Models Available

ChatBox40

The standard ChatBox40 is available in red or gray and includes both 10 and 40 location keyboards and keyguards. All ChatBox40 units have 30 minutes of digitized speech and have 10 overlays for programming messages. All messages in the ChatBox40 are spoken by pressing a single button.

ChatBox40 (Red) #2600RED\$895

ChatBox40-XT All of the features of the standard ChatBox40, with the addition of speech synthesis (male and female voices). ChatBox40-XT also includes auditory scanning and a serial output connector for keyboard emulation. The Carrying Case is included with the ChatBox40-XT. ChatBox40-XT (Red) #2601RED\$1195

ChatBox40-XT (Gray) #2601GRAY\$1195

S. ((1/2) B. DOCK P. " The ChatBox40® has helped Jessica communicate her needs. She **loves** it, adapted to it very easily, and it has helped eliminate some of her frustrations.

Jessica's Mother

ChatBox40Plus

ChatBox40Plus has all of the features of the ChatBox40, but adds the capability to store messages under a sequence of two icons, allowing instant access to many additional messages. The TalkToday vocabulary program, written by Gail VanTantenhove, is included in all Plus models.

ChatBox40-Plus (Red) #2615RED\$1095 ChatBox40-Plus (Gray) #2615GRAY\$1095

ChatBox40-XTPlus

The ChatBox40-XTPlus combines the speech synthesis and auditory scanning features of the ChatBox40-XT with the icon sequencing ability of the ChatBox40-Plus. The Carrying Case is included with the ChatBox40-XTPlus. ChatBox40-XTPlus (Red)

#2616RED\$1395 ChatBox40-XTPlus (Gray) #2616GRAY\$1395

Accessories

Carrying Case #2608\$65 Carrying case with shoulder strap (Carrying Case included with ChatBox40-XT and ChatBox40-XTPlus) Serial Cable #2609\$25 (Connects ChatBox40-XT to computer for keyboard emulation)









ChatBox® Simple, practical, effective

Are you searching for a sturdy, practical and cost-effective communication device? The Saltillo ChatBox offers digitized speech and customized options normally found in more expensive products.

Functional for children and adults, the ChatBox is a favorite of speech therapists. Easier to use and less expensive than other products with similar features, the ChatBox can be used with various accessories including the UniChat:16[™] vocabulary program.

Model CB-1

- > Volume Knob
- > On/off Switch
- > 16 active areas 1 " x 1 " each
- > 10 minutes of speech
- > No restrictions on message length
- > 4 levels
- > Levels can be changed from
- keyboard or back panel
- > Single key hits or sequencing (Minspeak[®])
- > Icon prediction
- > Rechargeable battery
- > Single switch scanning, Linear and row/column,
- > Latched, direct or step scanning > Scanning time variable:
- 0.5 5 seconds
- > Includes 4 overlays, customizing stickers, charger, manual
- > 12 month warranty
- > Size: 5.75" x 7.5"
- > Weight 1.5 lbs

Model CB-1

#2770\$495

Model CB-DX

- > All of the features of Model CB-1 (above) PLUS:
- > Auditory scanning with earphone jack and volume control
- > 4 remote switch inputs
- > Keyboard can be set up for either: 16 active areas 1 " x 1 " each or 4 active areas 2.25" x 2.25" each
- Model CB-DX

#2777\$645

Accessories for ChatBox ChatBox Carrying Cases

Carrying cases have a clear access

- window that allows the ChatBox to be
- used while inside the carrying case. #2775\$39 (includes shoulder strap)

#2776\$39 (includes waist belt)

UniChat:16[™] for ChatBox[®]

UniChat: 16 is a vocabulary program designed to be used by individuals learning to communicate on the ChatBox. Based upon Prentke Romich Company's Unity® program, UniChat: 16 is designed to give the system facilitator a starting point to organize vocabulary. It includes an instructional manual and 10 overlays to accommodate a variety of activities and settings.

UniChat:16

- #2772\$39 Set of 4 blank overlays (specify 4 or 16 locations)
- #2778\$2.50 Set of 20 overlays on 8.5 " x 11 " paper
- #2779\$7.50 (compatible with Mayer-Johnson Boardmaker[®] program)



ChatBox®

Model CB-1



Includes 4- and 16-location keyguards



VocaFlex[™] by Falck Igel AS Technology that makes magic

Simplicity and exceptional voice quality make the VocaFlex an outstanding product for use by young children and others in a group setting. The state-of-the-art voice quality is only surpassed by the unique and innovative technology behind it. Using bar coded overlays in conjunction with sturdy keyguards to form SmartGuides, the device requires virtually no training for you or the user. A canvas carrying case is available for easy transport of the VocaFlex.

The outstanding feature of the VocaFlex is the ability to automatically change the keyboard configuration and vocabulary by merely changing the overlay. The keyboards are attached by magnets so that they can be changed without assistance, giving the individual a greater sense of independence.

- > 3 position volume switch
- > Automatic on/off and a rechargeable Lithium-Ion battery
- > 9 FingerGuides included (3 each of 4, 8 and 16 locations)
- > 30 minutes or 1 hour of speech
- > Message lengths: 10 and 20 sec.
- > 51 Levels
- > Direct selection > Includes 90 overlays, vocabulary tracking chart and user's manual
- > 3 additional overlays for speaking the entire contents of an overlay, erasure of overlay content or erasure of all content of the device
- > 12 month warranty
- > Keyboard area: 3" x 3"
- > Size: 6.3" x 4" x 1.1"
- > Weight: 11 oz

VocaFlex with 30 minutes recording time #2881\$1695 VocaFlex with 1 hour recording time #2882.....\$1895



VocaFlex Accessories Canvas Carrying Case

A padded carrying case holds the VocaFlex and ten SmartGuides. VocaFlex can be used while in the carrying case. Shoulder strap and waist belt included. Canvas Caso



Callvas Case	
#2889	\$79
Additional SmartGuides	
4 Location (set of 5)	
#2886	\$135
8 Location (set of 5)	
#2887	\$135
16 Location (set of 5)	
#2888	\$135



Mounting Systems

FlexABLE by FRS Custom Solutions

FlexABLE mounts are a great solution for positioning ChatPC devices for easy access. These systems provide the flexibility and strength you need without the bulk and cost of mounting systems designed for much larger devices.

FlexABLE Stand and Adaptor

The FlexABLE Stand and Adaptor serves as a very flexible stand for a ChatPC-4 or ChatPC-Silk. This innovative new design provides infinite adjustment that allows the ChatPC to be positioned on a table or other flat surface for ease of access and glare reduction. The FlexABLE Stand and Adaptor can be used in either a portrait or landscape mode. The stand also serves as a convenient handle for transporting.

In addition to serving as a stand, the FlexABLE Stand and Adaptor can be attached to any of the mounting systems found on page 29.

FlexABLE Stand and Adaptor for ChatPC-Silk: #FLX-Silk\$179

FlexABLE Stand and Adaptor for

ChatPC-4: #FLX-4\$179





Adaptor for ChatPC-4

Mounting Systems for FlexABLE

The mounting systems below attach to a tube, table, wheelchair etc. They also attach to the FlexABLE Stand and Adaptors shown on page 28.

Tube Mount:

The tube mount attaches to a 7/8" tube. The rubber pads can be removed to fit a 1" tube. #FLX-TUBE\$129

Variable Friction Mount:

The variable friction mount can be manipulated into an infinite number of positions. One large, easy to grip knob tightens all of the joints, leaving the mount positioned securely. #FLX-VF\$370

Articulating Mount:

A lighter duty version of the Variable Friction Mount. The variable friction mount can be manipulated into an infinite number of positions. Each swivel joint is individually tightened.

#FLX-ART\$270



FlexABLE Stand and Adaptor for ChatPC-Silk attached to Articulating Mount









a 7/8" bar

Variable Friction Mount with FlexABLE Stand and Adaptor for ChatPC-4



Tube Mount with FlexABLE Stand and Adaptor for ChatPC-4

Note: A FlexABLE mounting solution consists of two components. Choose one of the mounting systems above. Also choose the proper FlexABLE stand and adaptor from page 28.





Access[™] Switches Control is only a touch away

Operate a communication aid, toy or environmental control with the new Saltillo Corporation Access Switches. These new switches are suitable for any device that uses one or more single input switches (standard 3.5mm plug). Access Switches are high quality products at a reasonable price and carry a limited lifetime warranty.

Designed to emit an audible "click" when activated, these switches require only a light touch. Constructed out of a high impact plastic, they are available in bright eye-catching colors. Measuring 2.5" in diameter, an individual can touch anywhere on the surface to activate it.

To provide a visible cue, images or symbols may be placed under the clear plastic cover which easily attaches to the top of the switch. Favorite photos can also be placed under the cover for a personal touch. The Access Switch design allows you to be both practical and creative!

Red Access Switch #2500\$29 Yellow Access Switch

#2501\$29 Blue Access Switch #2502\$29 Green Access Switch #2503\$29

Set of 4 Access Switches (1 of each color) #2505\$85



Access[™] Switches



Voice Amplifiers By Falck Igel AS The ultimate in voice quality

Place the Voice Amplifier in a pocket or around the neck – and those who speak in a whisper or in a quiet voice are ready to address the world.

Small, yet powerful, our top-of-the-line Voice Amplifier by Falck Vital offers fantastic voice reproduction at an affordable price.

The base of the amplifier is $2.5" \times 4.25"$ and weighs 4 oz. It can be worn around the neck using the adjustable "necklace", or it can be placed in one's shirt pocket. Membrane switches on the back of the amplifier are used to turn the unit on and off.

- > The Voice Amplifier is powered by rechargeable Lithium-Ion batteries.
- > The nominal usage for fully charged batteries is 12 hours. Recharging takes 4 hours.

Model 1003 Voice Amplifier includes microphone mounted on a headset. It also accepts the optional microphones. #3007\$975

Model 1001 Voice Amplifier has a



microphone mounted on a 7" swan-neck that can be adjusted for optimal placement.

#3006. ..\$975

Handheld Microphone for Model 1003 Amplifier

Clip-on Microphone for Model 1003 Amplifier

Clip-On Microphone

#3008. ..\$139

#3009.

.\$169

30 | www.saltillo.com



EchoVoice EV5™



EchoVoice EV5[™] By Hearing Products International Ltd. An all new compact voice amplification system that is sensitive enough to dramatically amplify voices.

EchoVoice EV5[™] is simple to use, comfortable, lightweight, discreet, and robust enough for classroom environments, sales meetings, group discussions or tour guides and many other areas where voice amplification is desirable.

The EV5™ is supplied with a lightweight headband type boom microphone and a transdermal (throat) microphone. Only a whisper is required to provide excellent amplification with the boom microphone. Alternatively the transdermal (throat) microphone requires positioning on the side of the throat, providing a more discreet option for the user.

A belt attachment on the casing and belt (provided) make the device easy to carry and simple to wear and remove.

Naphon Amplifier

The NAPHON voice amplifier is a cost-effective alternative for voice amplification. It features stable performance and good frequency response. Includes a belt and lightweight headset.

Device Specifications

- > Frequency Response: 50Hz-12KHz
- > Max. output power: 6W
- > Output impedance: 40hm
- > Operating Voltage: 6 AAA NI-MH rechargeable batteries
- > Charge time: 6-8 hours
- > Operating time: 4-6 hours
- continuous use
- Naphon Amplifier

#3095 \$95

. Сларноп Naphon Amplifier

Features

> Smart Carry Case

Controls

> Microphone and Master Volume

> Frequency Response 60 Hz – 12.5 kHz

> Total Harmonic Distortion < 0.5%

> Working Range Up to 50 yds

> Battery: Rechargeable Lithium

> Battery Run-Time: 5 to 8 hours

..\$595

#3090.....

> Intelligent Charge Circuit

> Excellent Battery Life

EchoVoice Specifications

> Dynamic Range 120 dB

> Power Output 10W

7.4V 1.2AmH

> Weight: 275 g

Accessories

On this page you will find various common accessories and replacement parts for Saltillo products. Some accessories are also listed on the individual product pages. Additional replacement parts and pricing for our discontinued products are available by calling 1-800-382-8622 or emailing service@saltillo.com.

Parts for NOVA chat 7:	Parts for ChatPC-Silk Product
NOVA chat 7 battery charger	ChatPC-Silk battery charger
#NOVA-BC\$69	#TD009
NOVA chat 7 carrying case	ChatPC-Silk USB cable (Used to conne
#NOVA-7-CC\$85	ChatPC-Silk to a desktop computer)
NOVA chat 7 shoulder strap	#WI015
#NOVA-7-SS\$25	ChatPC-Silk stylus
NOVA chat 7 waist belt	#HA023
#NOVA-7-WB \$25	ChatPC-Silk carrying case
NOVA chat 7 keyguard set (15/24/42)	#CC-Silk-2
#NOVA-7-KG Set\$60	ChatPC-Silk canvas carrying case
NOVA chat 7 keyguards (Individual)	#CC-CPC-Silk
#NOVA-7-(15)(24)(42)\$25	Parts for ChatPC-4 Products:
NOVA chat 7 wheelchair mounting plate	ChatPC-4 battery charger
#NOVA-7-WCMP \$59	#4006
NOVA chat 7 USB cable	ChatPC-4 USB cable (Used to connect
#NOVA-7-USB\$35	ChatPC-4 to a desktop computer)
NOVA chat 7 stylus	#WI010
#HA026\$10	ChatPC-4 stylus
Parts for NOVA chat 10:	#HA022
NOVA chat 10 battery charger	ChatPC-4 carrying case
#NOVA-BC\$69	#4005
NOVA chat 10 shoulder strap	
#NOVA-10-SS\$25	Parts for ChatBox-1 and
NOVA chat 10 keyguard (set of 7)	ChatBox-DX:
#NOVA-10-KG Set\$145	ChatBox battery charger
NOVA chat 10 keyguards (Individual)	#CBC-1
#NOVA-10-(specify size) \$35	Battery pack
NOVA chat 10 wheelchair mounting plate	(for serial numbers ending in A)
#NOVA-10-WCMP \$59	#SAL001BA
NOVA chat 10 USB cable	Battery pack
#NOVA-10-USB \$35	(for serial numbers ending in C or D)
NOVA chat 10 stylus	#SAL003BA
#HA026 \$10	Parts for ChatBox40:
	ChatBox40 battery charger
Parts for ALT-Chat:	#2604
ALT-Chat battery charger	ChatBox40 battery pack
#BC-PRC-12\$59	#BA007
ALT-Chat carrying case	Blank overlays (10 each of 10 and
#CC-SBL\$125	40 location)
ALT-Chat shoulder strap	#2605
#SS-ALT \$25	Printer overlays (10 each of 10 and
ALT-Chat stylus with lanyard	40 location)
#HA025\$15	#2606
Wheelchair mounting plate for ALT-Chat	#2000
#ALT-WCMP \$59	
36 Location PalmChat Keyguard for ALT-Chat	
#ALT-KG(36)\$125	
Custom keyguard for ALT-Chat	
(specify configuration)	
#ALT-KG(C) \$200	
Waist Belt for ALT-Chat	
#ALT-WB\$25	

Products:

....\$39 sed to connect computer) \$35 \$12 \$95 ng case \$65

Products:

Extended	Warranty
Extended	wanany

Parts for Falck Vital products:

#VF-1 \$20

#VF-BA \$80

#VA-1.....\$20

#EA006 \$20

VocaFlex battery charger

Voice amplifier battery charger

Parts for EchoVoice:

EchoVoice battery charger

VocaFlex battery

. \$39 Saltillo Corporation offers extended warranty on many of its products. The extended warranty has the same\$35 coverage as the original warranty. Up \$12 to four years of extended warranty may be purchased for a device. Also please\$65 note: extended warranty must be purchased before the original warranty, or previous extended warranty, expires. \$25 Extended Warranty for NOVA chat 10

	#EXW-NOVA-10	\$645/year
65	Extended Warranty for I	VOVA chat 7
	#EXW-NOVA 7	\$595/year
	Extended Warranty for A	ALT-Chat
65	#EXW-ALT	\$595/year
	Extended Warranty for (ChatPC-Silk
	#EXW-CPC-Silk	\$495/year
	Extended warranty for C	ChatPC-4
539	#EXW-CPC-4	\$495/year
	Extended warranty for Ch	atBox-1,ChatBox-DX
\$45	#EXW-CB	\$75/year
	Extended warranty for C	ChatBox-40
\$8	#EXW-CB40	\$150/year
	Extended warranty for V	/ocaFlex,
	#EXW-VF	\$150/year

.....\$12





How to Purchase Saltillo Products

Saltillo strives to make its products available to all who are in need of assistive technology. This is accomplished by working with reputable distributors on a national and international basis and partnering with Prentke Romich Company's Funding department to process purchases through third-party agencies such as Medicare, Medicaid and insurance companies.

For direct purchase information visit our online store www.silver-kite.com. There you can obtain pricing information and process an order with VISA, MC or American Express. We also accept prepayment by check or money order. Organizations may order with a Purchase Order. You may send your Purchase Orders to info@saltillo.com or fax to 330-674-6726.

For your convenience, Saltillo offers a monthly rental program. This program is ideal for determining if a device is appropriate for your communication needs or for installment payments for purchase of a device.

Find Your Regional Representative

Funding

Funding is available for your Saltillo device!

Helping you find the most appropriate technology for your client is very important to us. Saltillo representatives will assist you in the process from device trials to successful implementation strategies.

To find your regional representative go to www.saltillo.com.

E Store

Billy: H

Saltillo is pleased to introduce Silver-Kite, an assistive technology e-commerce website. Silver-Kite highlights Saltillo products as well as other products from leading manufacturers in the industry. At Silver-Kite you can download printer friendly product information sheets, get the latest pricing information or place your order.

Go to www.silver-kite.com

Get a Quote

Quotes are available from the Silver Kite online store. Just add products to your cart, and select "get quote" upon checkout.

To better serve our customers, Saltillo has partnered with the Funding Department at Prentke Romich Company to process requests for Medicaid, Medicare and private insurance. Prentke Romich Company's funding department is staffed with experienced professionals who help people with disabilities acquire the products and services that help them meet their needs to communicate. The PRC funding personnel are knowledgeable about any new requirements of the different agencies (Medicare, Medicaid or Insurance) and will provide the expertise needed to process a claim.

Please contact the PRC Funding Department directly by calling toll-free at 800-268-5224 or email at funding@prentrom.com . To access the PRC funding website go to http://funding.prentrom.com/



Order Form

Bill To:		Ship To: (if different)
Facility:		Facility:
Mailing Address:		Street Address:
City:		City:
State:	Zip:	State:
Phone: ()		Phone: ()
Fax: ()		Fax: ()

Description

Part Number

Shipping & Handling (USA) For orders totaling Add Up to \$35 \$7.00 \$35 to \$100 \$9.00 \$100 to \$1,000 \$16.00 Over \$1,000 no charge	Payment Method Check COD Purchase Order: Nu Visa MasterCard
Orders are shipped via ground. Call Saltillo Corporation for alternate	Expiration Date:
shipping or shipping outside of the USA.	Signature:

Mail or Fax to: Saltillo Corporation, 2143 TR112 Millersburg, OH 44654, USA (330) 674-6722 • 1 800 382-8622 • fax: (330) 674-6726

Please use product number and product description when placing orders.

Method and Terms of Payment: Prepayment by check or money order, COD, Visa/Mastercard. Organizations can purchase by a Purchase Order. Terms are net 30. Past due accounts will be assessed with a 1.5% per month (18% per year) penalty.

Warranty & Service: Warranty and service are provided by the manufacturers. Most products carry a one-year warranty. Saltillo Corporation will perform service on Saltillo, Unlimiter and Falck products. Service questions can be directed to service@saltillo.com.

Rental Units: Saltillo is able to provide a number of the items in this catalog as rental units. ChatBox rental is \$45/month. ChatBox40, VocaFlex, and voice amplifier rental is \$95/month. ChatPC-Silk and ChatPC-4 rentals are \$250/month. NOVA chat and ALT-Chat

new condition.

Medicaid and Medicare funding assistance available in most states.

J. (II UIIIEIE	11()		
/:			
Address:			
		Zip:	
: ()		
(
(Quantity	Price	Total
		Subtotal	
Imber			
	Sales Tax (0	hio Residents 6.5%)	
(3 digit code)	Shippir	ng and Handling	


rentals are \$350/month. Customers can rent for any number of months. All rental fees after the initial month will apply toward purchase of the unit. Contact Saltillo Corporation for a copy of the Rental Contract or download the Contract from www.saltillo.com.

TOTAL enclosed

Returns: Saltillo Corporation will accept return of any product within 30 days of the shipping date, provided the customer makes prior arrangements with Saltillo Corporation and receives a return authorization number. Materials must be returned prepaid and in

For Additional Information: Call us at 1-800-382-8622, email us at service@saltillo.com or visit our website at www.saltillo.com.

Price Changes: Prices and availability of products in this catalog may change without notice.





Saltillo

Saltillo Corporation 2143 TR112 Millersburg, OH 44654, USA

(330) 674-6722 **1 800 382-8622** fax: (330) 674-6726 email: info@saltillo.com

order online: www.silver-kite.com

In addition to the devices found in this catalog, products from various other manufacturers are available from Saltillo Corporation.

Visit our website at **www.saltillo.com** for more information.





0

STATISTICS.

10